

KO (Knee Orthosis)

General Description:

The term KO is an acronym that stands for Knee Orthosis and describes the part of the body that this device encompasses. This device extends from approximately mid calf to mid thigh level.

There are many categories of KOs: Metal designs and plastic and metal designs. The metal design consists of a metal structure shaped to the limb and upholstered with leather at the points where the device makes contact with a person's body. The plastic and metal design is the one most frequently encountered today and is usually a plastic device custom molded to the person's body with metal components in key structural areas only. The reason for the increased use of the plastic and metal design is that it is lighter in weight and is considered to be more cosmetic.

Application:

Since there are many different reasons to use a KO, there are many different designs. Each design has its own special features and its own specific way to be put on properly. The orthotist providing the device will instruct the patient on the proper way to put the KO on to make best use of the design features of the device. Also refer to the instruction pamphlet that is supplied with the device.

Cleaning and Maintenance:

The best way to clean a KO is to spray the inside with rubbing alcohol and wipe dry to remove body oils and residue. It can also be cleaned by wiping it out with a damp towel and anti-bacterial soap or anti-bacterial moist towelettes. Be sure to remove any excess soap as this can cause irritation to the skin.

Many KO's are designed to be used in or out of water. Check with the manufacturer's instructions to see if water is harmful to the device. Keep the KO away from excessive heat to prevent damage to the plastic.

Bannockburn 2101 Waukegan Road, Suite 104 Bannockburn, IL 60015 847-444-0690	Lincoln Park/Chicago 2551 N. Clark St., Suite 200 Chicago, IL 60614 773-472-3663	Rush/Chicago 1725 W. Harrison St., Suite 220 Chicago, IL 60612 312-942-2011	UIC/Chicago 1740 W. Taylor St. Room C100 Chicago, IL 60612 312-996-6450	Hickory Hills 8641 W. 95 th St. Hickory Hills, IL 60457 708-599-8336
Naperville 1551 Bond St., Suite 111 Naperville, IL 60563 630-637-4638	Oakbrook Terrace 1S376 Summit Ave., Ct. E Oakbrook Terrace, IL 60181 630-424-0392	Oak Park 401 Harrison Street Oak Park, IL 60304 708-383-2257	Schaumburg 1701 E. Woodfield Rd., Ste 555 Schaumburg, IL 60173 847-619-1701	Schererville 6629 W. Lincoln Hwy. (Rt. 30), Suite 1 PO Box 507 Schererville, IN 46375 219-864-9501

Tips and Problem Solving:

KO's can be worn directly against the skin for better suspension or over a cotton or neoprene sleeve for greater comfort.

Usually, the back straps just below the knee provide the most suspension on the calf. Be sure these straps are kept tight.

If there are any problems or questions, please contact the orthotist.



<p>Bannockburn 2101 Waukegan Road, Suite 104 Bannockburn, IL 60015 847-444-0690</p>	<p>Lincoln Park/Chicago 2551 N. Clark St., Suite 200 Chicago, IL 60614 773-472-3663</p>	<p>Rush/Chicago 1725 W. Harrison St., Suite 220 Chicago, IL 60612 312-942-2011</p>	<p>UIC/Chicago 1740 W. Taylor St. Room C100 Chicago, IL 60612 312-996-6450</p>	<p>Hickory Hills 8641 W. 95th St. Hickory Hills, IL 60457 708-599-8336</p>
<p>Naperville 1551 Bond St., Suite 111 Naperville, IL 60563 630-637-4638</p>	<p>Oakbrook Terrace 1S376 Summit Ave., Ct. E Oakbrook Terrace, IL 60181 630-424-0392</p>	<p>Oak Park 401 Harrison Street Oak Park, IL 60304 708-383-2257</p>	<p>Schaumburg 1701 E. Woodfield Rd., Ste 555 Schaumburg, IL 60173 847-619-1701</p>	<p>Schererville 6629 W. Lincoln Hwy. (Rt. 30), Suite 1 PO Box 507 Schererville, IN 46375 219-864-9501</p>
