



WHO (Wrist-Hand Orthosis)

General Description:

WHO is an acronym for Wrist-Hand Orthosis. A WHO is used for a variety of conditions. It can be used for positioning/stabilization, for range of motion assistance, or to stop a joint from flexing or extending excessively. It can also be used for fracture management. Materials can include canvas, metal, and low temperature thermoplastics. The physician determines use of the device. It may be worn all day or for nighttime use only.

Application:

- Loosen straps, place hand in device, putting thumb and fingers in their designated slots.
- Secure closures and adjust for comfort.

Cleaning and Maintenance:

The WHO can be cleaned by hand washing with mild soap and water. It must be air-dried.

Tips and Problem Solving:

If the wearer experiences pain, swelling, sensation changes, or unusual reactions, check the position and tightness of the device. If problems continue, consult the orthotist.



Bannockburn

2101 Waukegan Road,
Suite 104
Bannockburn, IL 60015
847-444-0690

Naperville

1551 Bond St., Suite 111
Naperville, IL 60563
630-637-4638

Lincoln Park/Chicago

2551 N. Clark St., Suite 200
Chicago, IL 60614
773-472-3663

Oakbrook Terrace

1S376 Summit Ave., Ct. E
Oakbrook Terrace, IL 60181
630-424-0392

Rush/Chicago

1725 W. Harrison St.,
Suite 220
Chicago, IL 60612
312-942-2011

Oak Park

401 Harrison Street
Oak Park, IL 60304
708-383-2257

UIC/Chicago

1740 W. Taylor St.
Room C100
Chicago, IL 60612
312-996-6450

Schaumburg

1701 E. Woodfield Rd.,
Ste 555
Schaumburg, IL 60173
847-619-1701

Hickory Hills

8641 W. 95th St.
Hickory Hills, IL 60457
708-599-8336

Schererville

6629 W. Lincoln Hwy.
(Rt. 30), Suite 1
PO Box 507
Schererville, IN 46375
219-864-9501