

## Humeral Fracture Brace

### General Description:

The Humeral Fracture Brace is a protective device worn on the upper arm to aid in healing and alignment of the injured bone.

### Application:

- Relax your arm and let it hang naturally. Sprinkle talcum powder on the upper arm and put on a clean stockinette (a special, open top sock). Smooth out any wrinkles.
- Apply the front section on the front of the upper arm, and then apply the back section. The back section must overlap the front section. If possible, have someone assist the person putting on the brace to get the best fit.
- Take the Velcro strap and pull it through the plastic loop, located in center. Next, pull the strap through the plastic loop on the back section, holding the back and front sections together snugly. Continue pulling the Velcro strap over the front until it fastens to itself. Take the remaining Velcro strap and pass it through the plastic loop on the back section and fasten it to itself. Equal pressure should be maintained on all straps. Straps may be adjusted to obtain a snug fit.

### Cleaning and Maintenance:

The best way to clean the brace is to spray the inside with rubbing alcohol and wipe dry to remove body oils and residue. It can also be cleaned by wiping it out with a damp towel and anti-bacterial soap, or anti-bacterial moist towelettes. Wash the cotton stockinette in cold water and air dry or machine dry on low.

Do not immerse the brace in water, as this will harm it.

Keep the brace away from excessive heat to prevent damage to the plastic.

---

<b>Bannockburn</b> 2101 Waukegan Road, Suite 104 Bannockburn, IL 60015 847-444-0690	<b>Lincoln Park/Chicago</b> 2551 N. Clark St., Suite 200 Chicago, IL 60614 773-472-3663	<b>Rush/Chicago</b> 1725 W. Harrison St., Suite 220 Chicago, IL 60612 312-942-2011	<b>UIC/Chicago</b> 1740 W. Taylor St. Room C100 Chicago, IL 60612 312-996-6450	<b>Hickory Hills</b> 8641 W. 95 <sup>th</sup> St. Hickory Hills, IL 60457 708-599-8336
<b>Naperville</b> 1551 Bond St., Suite 111 Naperville, IL 60563 630-637-4638	<b>Oakbrook Terrace</b> 1S376 Summit Ave., Ct. E Oakbrook Terrace, IL 60181 630-424-0392	<b>Oak Park</b> 401 Harrison Street Oak Park, IL 60304 708-383-2257	<b>Schaumburg</b> 1701 E. Woodfield Rd., Ste 555 Schaumburg, IL 60173 847-619-1701	<b>Schererville</b> 6629 W. Lincoln Hwy. (Rt. 30), Suite 1 PO Box 507 Schererville, IN 46375 219-864-9501

---

**Bathing:**

When the doctor has given permission to clean the arm or to change the stockinette, the arm may be cleaned as follows:

- Shower with the brace on.
- Follow the instructions above for removal of brace.
- Gently dry the arm and reapply the brace with a clean stockinette following the instructions above.

**Tips and Problem Solving:**

- With the braced arm, try to use the hand, wrist, elbow and shoulder as much as possible.
- The brace should be worn snugly, but not so tight as to be uncomfortable.
- Do not remove the brace unless instructed to do so by the doctor or orthotist.
- If there are any problems or questions regarding the brace, please contact the orthotist.




---

<p><b>Bannockburn</b> 2101 Waukegan Road, Suite 104 Bannockburn, IL 60015 847-444-0690</p> <p><b>Naperville</b> 1551 Bond St., Suite 111 Naperville, IL 60563 630-637-4638</p>	<p><b>Lincoln Park/Chicago</b> 2551 N. Clark St., Suite 200 Chicago, IL 60614 773-472-3663</p> <p><b>Oakbrook Terrace</b> 1S376 Summit Ave., Ct. E Oakbrook Terrace, IL 60181 630-424-0392</p>	<p><b>Rush/Chicago</b> 1725 W. Harrison St., Suite 220 Chicago, IL 60612 312-942-2011</p> <p><b>Oak Park</b> 401 Harrison Street Oak Park, IL 60304 708-383-2257</p>	<p><b>UIC/Chicago</b> 1740 W. Taylor St. Room C100 Chicago, IL 60612 312-996-6450</p> <p><b>Schaumburg</b> 1701 E. Woodfield Rd., Ste 555 Schaumburg, IL 60173 847-619-1701</p>	<p><b>Hickory Hills</b> 8641 W. 95<sup>th</sup> St. Hickory Hills, IL 60457 708-599-8336</p> <p><b>Schererville</b> 6629 W. Lincoln Hwy. (Rt. 30), Suite 1 PO Box 507 Schererville, IN 46375 219-864-9501</p>
--	--	--	---	---

---