

## **Below Knee Shrinker**

### **General Description:**

One of the initial steps in preparing your residual limb for a prosthesis is with the use of a shrinker. A shrinker is a tight, elastic stocking that equally compresses all parts of the residual limb. The purpose of the shrinker is to:

- 1) Reduce the amount of fluid or edema in the limb.
- 2) Keep the residual limb at a more consistent volume.
- 3) “Shape” the residual limb so that it is a more ideal shape and size to fit into a prosthesis.

### **Application:**

- 1) Using both hands, scrunch down the shrinker to the seam on the bottom and stretch the opening. The seam of the shrinker should be horizontal or from side to side.
- 2) Pull the shrinker up on the limb and pull up the sides as you travel up the limb.
- 3) The shrinker should completely cover the knee and go onto the thigh. The shrinker is very snug. If your limb is still healing, it may be difficult to put on initially. It may be easier to put on with assistance from another person (See Figure 1) or with the use of a donning can (See Figure 2).



Figure 1



Figure 2

**Arlington Heights**  
617 E. Golf Road, Suite 108  
Arlington Heights, IL 60005  
847-437-3929

**Bannockburn**  
2101 Waukegan Road,  
Suite 104  
Bannockburn, IL 60015  
847-444-0690

**Lincoln Park/Chicago**  
2551 N. Clark St., Suite 200  
Chicago, IL 60614  
773-472-3663

**Rush/Chicago**  
1725 W. Harrison St.  
Suite 220  
Chicago, IL 60612  
312-942-2011

**UIC/Chicago**  
1740 W. Taylor St.  
Room C100  
Chicago, IL 60612  
312-996-6450

**Hickory Hills**  
8641 W. 95<sup>th</sup> St.  
Hickory Hills, IL 60457  
708-599-8336

**Naperville**  
1551 Bond St., Suite 111  
Naperville, IL 60563  
630-637-4638

**Oakbrook Terrace**  
1S376 Summit Ave., Ct. E  
Oakbrook Terrace, IL 60181  
630-424-0392

**Oak Park**  
401 Harrison Street  
Oak Park, IL 60304  
708-383-2257

**Schererville**  
6629 W. Lincoln Hwy.  
(Rt. 30), Suite 1  
PO Box 507  
Schererville, IN 46375  
219-864-9501

- 4) Make sure the shrinker is pulled taut with no wrinkles. The shrinker must be pulled up all of the way at all times. This may require you to pull up the shrinker several times throughout a day to keep it snug. The shrinker may fall off while you are sleeping. If this happens, just remember to pull it back on as soon as you can.
- 5) Wear the shrinker at all times, day and night, except when you are bathing yourself or cleaning the shrinker. When you are cleaning the shrinker, a 2<sup>nd</sup> shrinker or an elastic bandage should be applied to the limb to prevent swelling.
- 6) Remember to keep your leg straight when sitting or laying down. Do not sit with your knee bent.

**Cleaning and Maintenance:**

Follow the manufacturer’s directions to clean the shrinker. Most recommend hand or machine wash on low settings with mild soap. Machine dry on low or hang to dry.

---

<p><b>Arlington Heights</b> 617 E. Golf Road, Suite 108 Arlington Heights, IL 60005 847-437-3929</p>	<p><b>Bannockburn</b> 2101 Waukegan Road, Suite 104 Bannockburn, IL 60015 847-444-0690</p>	<p><b>Lincoln Park/Chicago</b> 2551 N. Clark St., Suite 200 Chicago, IL 60614 773-472-3663</p>	<p><b>Rush/Chicago</b> 1725 W. Harrison St. Suite 220 Chicago, IL 60612 312-942-2011</p>	<p><b>UIC/Chicago</b> 1740 W. Taylor St. Room C100 Chicago, IL 60612 312-996-6450</p>
<p><b>Hickory Hills</b> 8641 W. 95<sup>th</sup> St. Hickory Hills, IL 60457 708-599-8336</p>	<p><b>Naperville</b> 1551 Bond St., Suite 111 Naperville, IL 60563 630-637-4638</p>	<p><b>Oakbrook Terrace</b> 1S376 Summit Ave., Ct. E Oakbrook Terrace, IL 60181 630-424-0392</p>	<p><b>Oak Park</b> 401 Harrison Street Oak Park, IL 60304 708-383-2257</p>	<p><b>Schererville</b> 6629 W. Lincoln Hwy. (Rt. 30), Suite 1 PO Box 507 Schererville, IN 46375 219-864-9501</p>

---